

Instruction for Use

【Product Name】

Troponin I (Electrochemiluminescence Immunoassay)

【Order Information】

REF NO.	Package Size
691056	50T
691057	2×50T
691011	100T
691012	2×100T

【Intended Use】

Immunoassay for in vitro quantitative determination of cardiac troponin I in human serum and plasma. This assay is intended to auxiliary diagnosis of myocardial infarction.

Cardiac troponin (Tn) was first reported by Cummins et al in 1987 when they diagnosed acute myocardial infarction (AMI)^{1,2}. Troponin (Tn) is the main structural contractile protein of myocardium, which is composed of three subunits: troponin C (TnC), troponin T (TnT) and troponin I (TnI)³. TnI can be divided into three subtypes: fast skeletal muscle subtype (fTnI), slow skeletal muscle subtype (sTnI) and myocardial subtype (cTnI). However, cTnI was the only one which could be detected in the myocardium when nine months after birth. cTnI has no subtype and is composed of 209 amino acids. More than 40% of the amino acid sequences are heterogenous to other TnI subtypes, and the amino terminal extends a 32-amino acid sequence, so its antigenicity is significantly different from skeletal muscle.

Clinical trials have shown that cTnI can be released into the blood system within hours of myocardial infarction (AMI) or ischemic injury^{4,5}. Elevated cTnI levels (above non-AMI sample values) can be detected in serum within 4-6 hours after the onset of chest pain, with the concentration peaking within 8-28 hours and remaining elevated for 3-10 days after the onset of AMI.

【Test Principle】

Sandwich principle. Total duration of assay: 18 minutes.

1st incubation: 100 µL of sample, biotinylated monoclonal cTnI-specific antibody and monoclonal cTnI-specific antibody labeled by ruthenium trisbipyridyl NHS ester form a sandwich compound by immunoreaction.

2nd incubation: After addition of streptavidin-coated micro magnetic beads, the compounds bind to the micro beads via interaction of biotin and streptavidin. Measurement: The reaction mixture is aspirated into the measuring cell where the micro beads and immuno-compounds are magnetically captured onto the surface of the electrode.

Unbound substances are then removed by Buffer. Application of a voltage to the electrode then induces emission of luminescence which is measured by a photomultiplier.

Results are determined via calibration and a master curve provided via the reagent barcode.

【Main Components】

The reagent pack consists of MB, RA, RB, and calibrators. Different lots cannot be used at the same time or mixed up together.

Components	Ingredients	Volume (50T)	Volume (2×50T)	Volume (100T)	Volume (2×100T)
(MB)	Streptavidin-coated micro magnetic beads 0.75mg/mL; 0.1 M PBS; 0.05% ProClin™ 300.	1×3.0 mL	2×3.0 mL	1×5.0 mL	2×5.0 mL
(RB)	Anti-cTnI-Ab~biotin 0.8 mg/L; 0.05 M MES buffer; 0.05% ProClin™ 300.	1×3.5 mL	2×3.5 mL	1×6.0 mL	2×6.0 mL
(RA)	Anti-cTnI-Ab~Ru(bpy) ₃ ²⁺ 0.63 mg/L; 0.05 M MES buffer; 0.05% ProClin™ 300.	1×3.5 mL	2×3.5 mL	1×6.0 mL	2×6.0 mL
Calibrators	Equine serum, 0.05% ProClin™ 300, lyophilized	Low: 1×2.0 mL High: 1×2.0 mL			

Traceability: This method can be traced back to NIST/SRM 2921.

Supporting instruments and materials required but not provided in this pack (provided by Lifotronic)

Additional materials for Automated ECL Immunoassay Analyzer eCL8000, eCL8000i, eCL8000p, eCL8000x:

- [REF] 679005, Auffer, 480 mL, or [REF] 679004, Auffer, 6×480 mL
- [REF] 679007, Buffer, 480 mL, or [REF] 679006, Buffer, 6×480 mL
- [REF] 679009, Concentrated Washing Buffer, 1 L, or [REF] 679008, Concentrated Washing Buffer, 6×1 L
- [REF] 679011, Assay Cup, 300 T, or [REF] 679012, Assay Cup, 3000 T
- [REF] 691027, Multi Control Cardiac Marker, Low value: 3×2.0 mL, High value: 3×2.0 mL, or [REF] 691028, Multi Control Cardiac Marker, Low value: 6×2.0 mL, High value: 6×2.0 mL

Additional materials for Automated ECL Immunoassay Analyzer eCL8600, eCL8600i, eCL8600p, eCL8800, eCL8800i, eCL8800p:

- [REF] 0320501804, Auffer, 2×1 L, or [REF] 0320501805, Auffer, 4×1 L
- [REF] 0320501904, Buffer, 2×1 L, or [REF] 0320501905, Buffer, 4×1 L
- [REF] 0320507502, PreClean, 2×800 mL, or [REF] 0320507503, PreClean, 4×800 mL
- [REF] 0320510001, Assay Cup/Assay Tip, 6×6×105, or [REF] 032120111,

Assay Cup/Assay Tip, 12×2×105

- [REF] 691027, Multi Control Cardiac Marker, Low value: 3×2.0 mL, High value: 3×2.0 mL, or [REF] 691028, Multi Control Cardiac Marker, Low value: 6×2.0 mL, High value: 6×2.0 mL

Additional materials for Automated ECL Immunoassay Analyzer eCL9000, eCL9600, eCL9900, eCL9900i:

- [REF] 0320501801, Auffer, 2×2 L, or [REF] 0320501806, Auffer, 4×2 L
- [REF] 0320501901, Buffer, 2×2 L, or [REF] 0320501906, Buffer, 4×2 L
- [REF] 0320507501, PreClean, 2×2 L, or [REF] 0320507504, PreClean, 4×2 L
- [REF] 0320510001, Assay Cup/Assay Tip, 6×6×105, or [REF] 0320510002, Assay Cup/Assay Tip, 1×6×105

- [REF] 691027, Multi Control Cardiac Marker, Low value: 3×2.0 mL, High value: 3×2.0 mL, or [REF] 691028, Multi Control Cardiac Marker, Low value: 6×2.0 mL, High value: 6×2.0 mL

【Storage and Shelf Life】

Store at 2°C-8°C.

Do not freeze.

Store the reagent kit upright in order to ensure complete availability of the microparticles during automatic mixing prior to use.

Damaged, expired or contaminated reagents should be discarded.

Stability of the reagents	
Unopened at 2°C-8°C	18 months
Store at 2°C-8°C after opening	8 weeks
Stored on the analyzers at 2°C-15°C	8 weeks

The calibrators are stable up to the stated expiration date.

Stability of the calibrators	
Unopened store at 2°C-8°C	18 months
Store at 2°C-8°C after reconstituted	8 hours
Store at -25°C~-15°C after reconstituted	3 months (freeze only once)

【Applicable Instrument】

Automated ECL Immunoassay Analyzers: eCL8000, eCL8000i, eCL8000p, eCL8000x, eCL8600, eCL8600i, eCL8600p, eCL8800, eCL8800i, eCL8800p, eCL9000, eCL9600, eCL9900, eCL9900i.

【Specimen collection, handling and storage】

Human serum and plasma added with heparin lithium, heparin sodium, EDTA-3K and sodium citrate anti-coagulants are recommended. Blood samples should be collected by standard operation of venous puncture. After complete coagulation, the tangible component should be removed by centrifugation. The sample should avoid bubbling during testing. Lipid layer floating on the upper of the sample should be removed. Samples with severe hemolysis are not recommended. Samples would better to be tested in 8 hours after collection, otherwise, should be stored 2-8°C for no more than 24 hours or -25~-15°C, 3 months. Freezing and thawing cycle is permitted only once.

【Assay Procedure】

Testing procedures and precautions

Previous to operation, one should read the operational manual carefully to obtain system operation procedure, sample processing, security precaution, maintenance and other related information.

Set up Troponin I (cTnI) test according to the operational manual.

Put the cTnI rack pack into the correct analyzer slot to automatically suspend magnetic beads at least 30 min before testing.

Suggested testing environment is 10-30°C; Relative humidity less than 80%.

Handling of lyophilized calibrators:

Carefully dissolve the contents of one bottle by adding exactly 2.0 mL of distilled or deionized water and allow to stand closed for 10-15 minutes to reconstitute. Mix carefully, avoiding foam formation. Transfer aliquots of the reconstituted calibrators to the additional bottles. After bottles is marked and dispensed. If you do not test immediately, you need to store the aliquots immediately at -25~-15°C.

Calibration

Calibration should be performed using lot-matching reagent and calibrators.

Before calibration, the reagent and main curve information should be imported into the analyzer via radio frequency identification (RFID) reagent card (refer to instrument operational manual). The analyzer adjusts the main curve to produce working curve according to the calibrator results, and identifies validity of the working curve automatically.

Recalibration is recommended when:

- (1) The reagents of different lot are used;
- (2) The same lot of reagents were used on the analyzer beyond 4 weeks;
- (3) Quality control misses the target.
- (4) The Buffer of different lots are used.

Quality control

For quality control, use Multi Control Cardiac Marker.

In addition, other suitable control material can be used.

Control materials for the various concentration ranges should be run individually at least once every 24 hours when the test is in use, once per reagent kit, and

following each calibration.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

If necessary, repeat the measurement of the samples concerned.

Follow the applicable government regulations and local guidelines for quality control.

Please make sure that the correct values are used.

Calculation

The system software automatically calculates the analyte concentration using particular algorithm, and the result unit is ng/mL, µg/L or pg/mL:

1 ng/mL = 1 µg/L; 1 ng/mL = 1000 pg/mL

Specimen dilution

Samples with cTnI concentrations above the measuring range can be diluted with cTnI-negative sample diluent. The recommended dilution ratio is 1:10. The concentration of the diluted sample must be > 3 ng/mL. After manual dilution, multiply the result by the dilution factor.

【Biological reference interval】

According to 312 cases of serum samples of healthy adult without heart disease human in Guangdong hospitals, P. R. China, the 99 percentile draws a reference range as follow:

Reference range	cases	99 percentile (ng/mL)
10~85 years old	312	0.028

Acute myocardial infarction (AMI) cut-off value

By analyzing 203 patients with chest pain, according to the world health organization (WHO) AMI diagnostic criteria, of which 73 for AMI, we set up cut-off value of AMI 0.1 ng/mL using receiver-operating characteristic (ROC) curve method.

As the differences in geography, race, gender and age, each laboratory should investigate transferability of the reference range to its local population and if necessary determine its own reference ranges.

【Result Interpretation】

In interpreting the results, the patient's overall clinical situation should be referred to, including symptoms, medical history and other relevant data and information.

【Limitations】

Test results are used only for clinical reference and cannot be used as the basis for diagnosis or rule-out of diseases alone.

The measuring range of the kit is 0.006~100 ng/mL. Values below lower detection limit are reported as < 0.006 ng/mL. Values above the measuring range, the sample need be diluted and the results will be reported as the measured value multiplied by the dilution factor.

When samples contain the triglyceride concentration of 1000 mg/dL or less, bilirubin concentration of 20 mg/dL or less, hemoglobin concentration of 500 mg/dL or less, total protein concentration of 8 g/dL or less, rheumatoid factor concentration of 1500 IU/mL or less, biotin concentration of 20 ng/mL or less, HAMA concentration at 100 ng/mL or less, the interference bias of determination results deviation within ± 15%.

There is no high-dose hook effect at cTnI concentration 2000 ng/mL.

【Analytical Performance】

Lower limit of measurement

The lower detection limit is 0.006 ng/mL (repeatability study, n = 20).

Accuracy

When cTnI of known concentration was added to the matrix without the tested substance, the recovery rate should be within (85%~115%).

Linearity

The correlation coefficient (r) was not less than 0.9900 in the interval of 0.006~100 ng/mL.

Within-run precision (repeatability)

The coefficient of variation (CV) ≤ 6%.

Between-lot precision

The coefficient of variation (CV) ≤ 10%.

Analytical specificity

Test the following analog-additive blank samples, the cross reaction rate as below.

Analog	Concentration	Cross reaction rate (%)
Sk-TnI	1000 ng/mL	0.001%
cTnC	1000 ng/mL	0.001%
cTnT	1000 ng/mL	0.002%

Method Comparison

A comparison of the Lifotronic cTnI assay (y) with the Abbot high-sensitivity troponin-I assay testing kit (x) using clinical samples gave the following correlations:

Number of serum samples measured:110

Linear regression:

y=0.9944x -0.0620

r = 0.996

The sample concentrations were between 0.013 and 105.261 ng/mL.

【Biosafety】

The HBsAg, HCV antibody and HEV antibody test results of calibrators in the kit are all negative. Nevertheless, any test could not exclude absolutely possibility of contagion. Thus, essential and regulation-compliant precautions should be adopted when operating calibrators in the kit.

【Precaution and Warning】

The kit is only used for in vitro diagnostics.

When using the kit, it's necessary to comply with regulations in the laboratory.

All reagents and samples including specimen, calibrators and control materials, should avoid foaming before and during test.

Test results of the kit are for clinical reference only, and clinical evaluation of patients should be combined with their symptoms and signs, medical history, other laboratory examination results and treatment responses.

Due to factors such as methodology or antibody specificity, testing identical samples with reagents from different manufacturers may get different results, and the results from different kits should not compare directly, lest cause the wrong medical explanation. It's suggested that laboratorians should point out characteristics of the used reagent in the test report to the clinician. In serial monitoring, if the reagent type is changed, a continuous parallel test should be performed and comparison of the results to the former to determine new baseline value.

This product contains animal-sourced materials and may have potential biological risk. All samples and test wastes should be treated as the source of infection, and all wastes must be disposed according to local regulations.

The preservative ProClin™300 which contain Mixture of 5-Chloro-2-methyl-4-isothiazolin-3-one and 2-Methyl-2H-isothiazol-3-one (3:1): May cause an allergic skin reaction. thus, proper personal protection should be adopt while handling the reagent: avoid breathing dust/fume/gas/mist/vapours/spray, wear protective gloves and contaminated work clothing should not be allowed out of the workplace. If skin irritation or rash occurs: get medical advice/attention, take off contaminated clothing and wash it before reuse, and dispose of contents/container to an approved waste disposal plant.

【Symbol Explanation】

Symbol	Title	Symbol	Title
	Manufacturer		Consult instructions for use
	Use-by date		In vitro diagnostic medical device
	Lot number		Contains sufficient for <n> tests
	Serial number		Biological risks
	Temperature limit		This way up
	Authorized representative in the European Community		Indicates this device is in compliance with Europe Directive.
	Catalogue number		Material code
	Warning		Date of manufacture
	Unique device identifier		

【Literature References】

- Cummins P, Perry V. Troponin I from human skeletal and cardiac muscles. Biochem J 1978;171:251-259.
- Wilkinson JM, Grand RJA. Comparison of amino acid sequence of troponin I from different striated muscles. Nature 1978;271:31-35.
- Mair J, Morandell D, Genser N, et al. Equivalent early sensitivities of myoglobin, creatine kinase MB mass, creatine kinase isoform ratios, and cardiac troponins I and T for acute myocardial infarction. Clin Chem 1995;41:1266-1272.
- Tanasijevic MJ, Cannon CP, Antman EM. The role of cardiac troponin I (cTnI) in risk stratification of patients with unstable coronary artery disease. Clin Cardiol 1999;22:13-16.
- Jeffrey L. Anderson, Cynthia D. Adams, Elliott M. Antman, et al. 2011 ACCF/AHA Focused Update of the Guidelines for the Management of Patients With Unstable Angina/Non-ST-Elevation Myocardial Infarction (Updating the 2007 Guideline): A Report of the American College of Cardiology Foundation/American Heart Association Task Force on Practice Guidelines. Circulation 2011.

【Manufacturer】

Shenzhen Lifotronic Technology Co., Ltd.

Lifotronic Tower, No.8, Qiuzhi East Road, Guancheng Community, Guanhu Street, Longhua, 518110 Shenzhen, PEOPLE'S REPUBLIC OF CHINA

E-mail: inter-service@lifotronic.com

【European Representative】

Shanghai International Holding Corp. GmbH (Europe)

Eiffelstrasse 80, 20537 Hamburg, Germany

Tel. +49-40-2513175 Fax. +49-40-255726

【Version and Revision】

Version: 2.0

Issue Date:2025-03-04