

Instruction For Use

【Product Name】

Free Thyroxine Kit (Electrochemiluminescence Immunoassay)

【Package】

REF NO.	Package Size
690052	50T
690053	2×50T
690054	100T
690055	2×100T

【Intended Use】

Immunoassay for in vitro quantitative determination of free thyroxine in human serum and plasma.

Thyroxine (T4) is the main thyroid hormone secreted into the bloodstream by the thyroid gland¹. Together with triiodothyronine (T3) it plays a vital role in regulating the body's metabolic rate, and has influence on the cardiovascular system, growth and bone metabolism, and as well as normal development of gonadal functions and nervous system².

The determination of free T4 has an advantage of independency of changes in the concentrations and binding properties of thyroxine-binding proteins; additional determination of a binding parameter (T-uptake, TBG) is therefore unnecessary³. Therefore, free T4 is a useful tool in clinical routine diagnostics for the assessment of the thyroid status. It should be measured together with TSH if suspected thyroid disorders and is also suitable for monitoring thyrostatic therapy⁴.

【Test Principle】

Competition principle. Total duration of assay: 18 minutes

15µL of sample and T4-specific antibody labeled with a ruthenium complex react to form a hapten-antibody complex. After addition of biotinylated T4 and streptavidin-coated microparticles, biotinylated T4 derivative competitively occupied the binding-free antibody and bound to the solid phase via interaction of biotin and streptavidin. The reaction mixture is aspirated into the measuring cell where the microparticles are magnetically captured onto the surface of the electrode. Unbound substances are then removed with Buffer. Application of a voltage to the electrode then induces chemiluminescent emission which is measured by a photomultiplier. Results are determined via calibration and a master curve provided via the reagent barcode.

【Main Component】

The reagent pack consists of MB, RA, RB, calibrators and control materials (Optional). Different batch of reagent components should not be used interchangeably.

Components	Ingredients	Volume (50T)	Volume (2×50T)	Volume (100T)	Volume (2×100T)
Magnetic Beads (MB)	Streptavidin-coated microparticles, 0.75mg/mL; 0.1M PBS; preservative	1×1.8 mL	2×1.8 mL	1×3.5 mL	2×3.5 mL
Reagent B (RB)	T4-Biotin, 1 ng/mL; 0.1 M PBS; preservative	1×3.5 mL	2×3.5 mL	1×7.0 mL	2×7.0 mL
Reagent A (RA)	T4-Ab-Ru(bpy) ₃ ²⁺ , 200 ng/mL; 0.1 M PBS; preservative	1×3.5 mL	2×3.5 mL	1×7.0 mL	2×7.0 mL
Calibrators	0.1 M PBS, preservative	High: 1×1.0 mL Low: 1×1.0 mL			
Control Materials (Optional)	0.1 M PBS, preservative	High: 1×1.0 mL Low: 1×1.0 mL			

Traceability: This method can be traced back to Roche Elecsys FT4 III.

See the control material value sheet for the target values and ranges of control materials.

Supporting instruments and materials required but not provided in this pack (provided by Lifotronic)

Additional materials for Automated ECL Immunoassay Analyzer eCL8000, eCL8000i, eCL8000p, eCL8000x:

- [REF](#) 679005, Auffer, 480 mL, or [REF](#) 679004, Auffer, 6×480 mL
- [REF](#) 679007, Buffer, 480 mL, or [REF](#) 679006, Buffer, 6×480 mL
- [REF](#) 679009, Concentrated Washing Buffer, 1 L, or [REF](#) 679008, Concentrated Washing Buffer, 6×1 L

- [REF](#) 679011, Assay Cup, 300 T, or [REF](#) 679012, Assay Cup, 3000 T

Additional materials for Automated ECL Immunoassay Analyzer eCL8600, eCL8600i, eCL8600p, eCL8800, eCL8800i, eCL8800p:

- [REF](#) 0320501804, Auffer, 2×1 L, or [REF](#) 0320501805, Auffer, 4×1 L
- [REF](#) 0320501904, Buffer, 2×1 L, or [REF](#) 0320501905, Buffer, 4×1 L
- [REF](#) 0320507502, PreClean, 2×800 mL, or [REF](#) 0320507503, PreClean, 4×800 mL
- [REF](#) 0320510001, Assay Cup/Assay Tip, 6×6×105, or [REF](#) 032120111, Assay Cup/Assay Tip, 12×2×105

eCL9600, eCL9900, eCL9900i:

- [REF](#) 0320501801, Auffer, 2×2 L, or [REF](#) 0320501806, Auffer, 4×2 L
- [REF](#) 0320501901, Buffer, 2×2 L, or [REF](#) 0320501906, Buffer, 4×2 L
- [REF](#) 0320507501, PreClean, 2×2 L, or [REF](#) 0320507504, PreClean, 4×2 L
- [REF](#) 0320510001, Assay Cup/Assay Tip, 6×6×105, or [REF](#) 0320510002, Assay Cup/Assay Tip, 1×6×105

【Storage and Shelf Life】

Store at 2°C-8°C.

Do not freeze.

Store the reagent kit upright in order to ensure complete availability of the microparticles during automatic mixing prior to use.

Stability of the reagents	
Unopened store at 2°C~8°C	18 months
Store at 2°C~8°C after opening	8 weeks
on the analyzers at 2°C~15°C	8 weeks

The calibrators and control materials are stable up to the stated expiration date.

Stability of the calibrators and control materials	
Unopened store at 2°C~8°C	18 months
Store at 2°C~8°C after opening	8 weeks

The expiration date is labeled on the box, rackpack and bottles.

Damaged, expired or contaminated reagents should be discarded.

【Applicable Instrument】

Automated ECL Immunoassay Analyzers: eCL8000, eCL8000i, eCL8000p, eCL8000x, eCL8600, eCL8600i, eCL8600p, eCL8800, eCL8800i, eCL8800p, eCL9000, eCL9600, eCL9900, eCL9900i.

【Specimen collection, handling and storage】

Human serum and plasma added with heparin lithium, heparin sodium and EDTA-K2 anti-coagulants are recommended. Blood samples should be collected by standard operation of venous puncture; after complete coagulation, the tangible component should be removed by centrifugation. The sample should avoid bubbling during testing. Lipid layer floating on the upper of the sample should be removed. Samples with severe hemolysis are not recommended. It is recommended to complete the test in time after sample collection. Samples can stable for 8 hours at 18~28°C, 7 days at 2~8°C, 28 days at -25°C~-15°C, Freeze only once.

【Assay Procedure】

Testing procedures and precautions

Previous to operation, one should read the operational manual carefully to obtain system operation procedure, sample processing, security precaution, maintenance and other related information.

Set up Free Thyroxine (FT4) test according to the operational manual.

Put the FT4 rackpack into the correct analyzer slot to automatically suspend magnetic beads at least 30 min before testing.

Calibration

Calibration should be performed using lot-matching reagent and calibrators.

Before calibration, the reagent and main curve information should be imported into the analyzer via radio frequency identification (RFID) reagent card. The analyzer adjusts the main curve to produce working curve according to the calibrator results, and identifies validity of the working curve automatically.

Recalibration is recommended when:

- (1) the lots of reagents are changed;
- (2) the same lot of reagents were used on the analyzer beyond 28 days;
- (3) quality control misses the target;
- (4) the lots of buffer are changed.

Quality control

For quality control, use FT4 control materials.

In addition, other suitable control material can be used.

Control materials for the various concentration ranges should be run individually at least once every 24 hours when the test is in use, once per reagent kit, and following each calibration.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

If necessary, repeat the measurement of the samples concerned.

Follow the applicable government regulations and local guidelines for quality control.

Please make sure that the correct values are used.

Calculation

The system software automatically calculates the analyte concentration using particular algorithm, and the result unit is pmol/L, ng/L or ng/dL: pmol/L × 0.077688= ng/dL;

ng/dL×12.872=pmol/L;

pmol/L×0.77688=ng/L.

Specimen dilution

Samples for FT4 determination can not be diluted, as T4 in the blood is present in free and protein-bound forms in equilibrium. Dilution will change free thyroxine

concentration by breaking the equilibrium.

【Biological reference interval】

According to serum test results of 294 seemingly healthy human (male 150, female 144), the percentile 2.5 to 97.5 gives reference range of 10.51~22.7 pmol/L.

Due to differences in the region, race, gender and age, laboratories are recommended to set up their own reference ranges.

【Result Interpretation】

For explaining the result, the patient's overall clinical manifestation should be referred to, including symptoms, medical history and other relevant data and information.

【Limitations】

Test results are used only for clinical reference and cannot be used as the basis for diagnosis or rule-out of diseases alone.

Deviation of results falls within $\pm 10\%$ with interference by triglyceride ≤ 1400 mg/dL, bilirubin ≤ 10 mg/dL, hemoglobin ≤ 500 mg/dL, total protein ≤ 1.15 g/dL, rheumatoid factor ≤ 800 IU/mL or biotin ≤ 10 ng/mL.

【Measuring Range】

The measuring range of the kit is 0.30~100 pmol/L, which is determined by lower limit of detection and upper limit of master calibration curve. Values below the lower detection limit are reported as < 0.3 pmol/L, values above the measuring range are reported as >100 pmol/L.

【Analytical Performance】

Lower limit of measurement

Analytical sensitivity is 0.30 pmol/L (repeatability study, $n = 20$).

Accuracy

The relative bias of measuring trueness control materials within $\pm 15\%$.

Linearity

The correlation coefficient (r) is not less than 0.99 in the interval of 0.30-100 pmol/L.

Within-run precision

The coefficient of variation (CV) $\leq 7.5\%$.

Analytical specificity

The following cross-reactivities were found.

Cross-reactant	Concentration	Cross-reactivity
L-T3	50,000 ng/dL	0.04100%
D-T3	50,000 ng/dL	0.04800%
rT3	950 ng/dL	0.00100%
3-iodo-L-tyrosine	10,000,000 ng/dL	0.00010%
3,5-diiodo-L-tyrosine	10,000,000 ng/dL	0.00100%
3,3',5,5'-tetraiodothyroacetic acid	100,000 ng/dL	0.00003%

Method comparison

A comparison of the Lifotronic FT4 assay (y) with the Roche Elecsys FT4

Method (x) using clinical samples gave the following correlations:

Number of serum samples measured: 233

Linear regression:

$$y=1.1617x-0.8725$$

$$r=0.975$$

The sample concentrations were between approx 2.43 and 94.34 pmol/L.

【Precaution and Warning】

The kit is only used for in vitro diagnostics.

When using the kit, it's necessary to comply with regulations in the laboratory.

All reagents and samples including specimen, calibrators and control materials, should avoid foaming before and during test.

Test results of the kit are for clinical reference only, and clinical evaluation of patients should be combined with their symptoms and signs, medical history, other laboratory examination results and treatment responses.

Due to factors such as methodology or antibody specificity, testing identical samples with reagents from different manufacturers may get different results, and the results from different kits should not compare directly, lest cause the wrong medical explanation; It's suggested that laboratorians should point out the characteristics of used reagent in the test report to the clinician. In serial monitoring, if the reagent type is changed, a continuous parallel test should be performed and comparison of the results to the former to determine new baseline value.

This product contains animal-sourced materials and may have potential biological risk. All samples and test wastes should be treated as the source of infection, and all wastes must be disposed according to local regulations.

Any serious incident that would be occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

The preservative ProClin™300 which contain Mixture of 5-Chloro-2-methyl-4-isothiazolin-3-one and 2-Methyl-2H-isothiazol-3-one (3:1): May cause an allergic skin reaction. thus, proper personal protection should be

adopt while handling the reagent. avoid breathing dust/fume/gas/mist/vapours/spray, wear protective gloves and contaminated work clothing should not be allowed out of the workplace. If skin irritation or rash occurs: get medical advice/attention, take off contaminated clothing and wash it before reuse, and dispose of contents/container to an approved waste disposal plant.

【Reference】

1.Pierce JG. Endocrinology 1971, 89:1331-44.

2.Morgan F and Canfield RC. Endocrinology 1971, 88(4): 1045-53.

3.Hertz R, Odell WD, and Ross GT. Ann Intern Med 1966, 65(4):800-20.

4.Fisher DA. Clin Chem 1996, 42(1):135-9.

【Symbols Explanation】

Symbol	Title of Symbol	Symbol	Title of Symbol
	Manufacturer		Consult instructions for use
	Use-by date		In vitro diagnostic medical device
	Batch code		Sufficient for $<n>$ tests
	Serial number		Biological risks
	Temperature limit		This way up
	Authorized representative in the European Community		Indicates this device is in compliance with Europe Directive.
	Catalogue number		Material code
	Warning		Date of manufacture
	Unique device identifier		

【Manufacturer】

Shenzhen Lifotronic Technology Co., Ltd.

Lifotronic Tower, No. 8, Qiuzhi East Road, Guancheng Community, Guanhu Street,

Longhua, 518110 Shenzhen, PEOPLE'S REPUBLIC OF CHINA

E-mail: inter-service@lifotronic.com

【European Representative】

Shanghai International Holding Corp. GmbH (Europe)

Eiffelstrasse 80, 20537 Hamburg, Germany

Tel. +49-40-2513175 Fax. +49-40-255726

【Version and Revision】

Version: 2.0

Issue Date: 2025-03-04